

Purification of RNA from Crude NA Extract

The following protocol is intended for the purification of DNA-free RNA from a crude extract of nucleic acids.

Reagents and kits

1. DNase I (preferably Ambion's TURBO DNase).
2. Qiagen's RNeasy MinElute Cleanup.
3. Optional: RNase inhibitor (Invitrogen's RNaseOUT™, Promega's RNasin™, or Fermentas' Ribolock™).
4. Optional: Ambion's RNA Storage Solution.

Before you begin

1. Optional: add 10 µL β-ME or 20 µL 2M DTT per 1 mL buffer RLT. Store up to 1 month at RT.
2. For each sample prepare in a rack: 1 MinElute spin column, 2 collection tubes, 1 1.5 mL tube and mark them.
3. Dilute buffer RPE with 44 mL of 96-100% EtOH.

Procedure

1. a. If using Qiagen's DNase I or Ambion's TURBO DNase prepare the following mixture in a 1.5 mL tube:

NA crude extract	10-85.5 µL
10X Buffer	10 µL
RNase inhibitor (40 U/µL)*	2 µL
0.1 M DTT*	1 µL
DNase I	1 µL per 2 µg of DNA or 2.5 µL for Qiagen's DNase I
RNase free water	To 100 µL

b. If using Fermentas' DNase I prepare the following mixture in a 1.5 mL tube:

NA crude extract	10-86 µL (or up to 1 µg per unit)
10X reaction buffer without MgCl ₂	10 µL
MnCl ₂ (100 mM)	1 µL
RNase inhibitor (40 U/µL)*	2 µL
0.1 M DTT*	1 µL
DNase I (1 U/µL)	1 µL per 1 µg of DNA
RNase free water	To 100 µL

* optional

2. Incubate at **37° C** (**room temperature** for Qiagen's DNaseI) for **30 min**.
3. Add **350 µL** buffer **RLT** (per 100 µL digested sample) and mix well.
4. Add **250 µL 96-100% EtOH** (per 100 µL digested sample) and mix well by pipetting.
5. Transfer immediately up to **700 µL** of sample to the spin column, close the lid, and centrifuge for **15 s at $\geq 8000 \times g$** (≥ 10000 rpm).
6. If sample volume is larger than **700 µL** discard the flow through after centrifugation and repeat the process until all the sample has been passed through the column.
7. Place the spin column in a new collection tube and add **500 µL** of buffer **RPE**, close the lid, and centrifuge for **15 s at $\geq 8000 \times g$** (≥ 10000 rpm).
8. Discard the flow through and place spin column in the same collection tube.
9. Add **500 µL** of **80% EtOH**, close the lid, and centrifuge for **2 min at $\geq 8000 \times g$** (≥ 10000 rpm).
10. Discard the flow through and place spin column in a new collection tube.
11. Open the lid, and centrifuge for **5 min at full speed**.
12. Place the spin column in a **1.5 mL** tube, add **20 µL** (10-30 µL) of **RNase free water** or **RNase Storage Solution**, and centrifuge for **1 min at full speed**.