

## Purification of DNA from Crude NA Extract

The following protocol is intended for the purification of RNA-free DNA from a crude extract of nucleic acids.

### Reagents and kits

1. RNase ONE™ Ribonuclease (Promega).
2. Wizard® DNA Clean-Up System (Promega).
3. Isopropanol.

### Before you begin

1. For each sample prepare in a rack: 1 – MinElute spin column, 2 – collection tubes, 1 – 1.5 ml tube and mark them.

### Procedure

1. Incubate **10-100 µl** of crude NA extract (up to 10 µg of total RNA at **85 °C** for **5 min** to denature the RNAs).
2. Prepare the following mixture in a 1.5 ml tube:

NA crude extract	50-89 µl (up to 20 µg)*
RNase digestion buffer	10 µl
RNase ONE™ Ribonuclease	1 µl
RNase free water	To 100 µl

\* If sample volume is less than 50 µl, bring to 50 µl with DDW.

3. Incubate the samples for **60 min** at **37 °C**.
4. Attach one Syringe Barrel to a Minicolumn, and insert the tip to the vacuum manifold.
5. Shake well the Wizard® DNA Clean-Up Resin and add **1 ml** to the sample tube and mix by inversion.
6. Pipette the resin/ sample mix into the Syringe Barrel.
7. Apply **vacuum** until solution is completely drawn.
8. Add **2 ml 80% isopropanol**, and re-apply **vacuum** until solution is completely drawn.

9. Continue applying the vacuum for extra **30 s** after the solution has been drawn.
10. Remove the Syringe Barrel and transfer the minicolumn to a 1.5 ml tube, and centrifuge for **2 min** at **maximum speed**.
11. Transfer the minicolumn to a new 1.5 ml tube, apply **50 µl** of **70 °C water** or **low TE buffer** and incubate **1 min**.
12. Centrifuge the minicolumn at **maximum speed** for **30 s**.
13. Run **2-5 µl** on agarose gel to validate removal of RNA.
14. Quantify the DNA using a spectrophotometer or PicoGreen.