16S Methanosarcinales 812F-860P-1159R

QPCR mixture and conditions

Name	Sequence	Target	Reference
MSL 812F	GTA AAC GAT RYT CGC TAG GT	812-831	Yu et al. B&B 2004
MSL 860P	AGG GAA GCC GTG AAG CGA RCC	860-880	u
MSL 1159R	GGT CCC CAC AGW GTA CC	1143-1159	u

Reagent	Final conc.	1 tube	Plate
		(25 μL)	(25 μL × 100)
PCR H ₂ O		1.6	160
ReadyMix (Sigma)	1X	12.5	1250
MgCl ₂ (25 mM)	4 mM	2.5*	250*
BSA (50 μg/μL)	0.8 μg/μL	0.4	40
812F (10 μM)	0.5 μΜ	1.25	125
1159R (10 μM)	0.5 μΜ	1.25	125
860P [†] (10 μM)	0.2 μΜ	0.5	50
Template		5	5 × 100

^{*}Buffer contains 1.5mM of MgCl₂ at final conc.

Programme:

1. 95 °C -5'

X40 {

2. 95 °C - 30"

3. $60 \,^{\circ}\text{C} - 60'' - \text{take snapshot}$

[†] Probe must be dual-labelled either with 5'-6-FAM, 3'-BHQ1 or any other valid combination.